



800-640-7780

Sensitivity Result Report

Patient Name _____

Date _____

Organism ID _____

Specimen Site _____

Zone Size Interpretation

Antibiotic		Resistant	Intermediate	Susepible	Results	circle one
Amikacin	AN	<= 14	15 - 16	>= 17		R S
Amoxicillin&Clavulanate =Augmentin	AUG	<= 13	14 - 17	>= 18		R S
Ampicillin	AM	<= 14	15 - 17	>= 18		R S
Azithromycin	ATH	<= 13	14 - 17	>= 18		R S
Cefazolin	CZ	<= 14	15 - 17	>= 18		R S
Cefovecin = Covenia®	CEF	<= 20	21 - 23	>= 24		R S
Cefoxitin	FOX	<= 14	15 - 17	>= 18		R S
Cefpodoxime	CPD	<= 17	18 - 20	>= 21		R S
Ceftazidime & Avibactam	CAZ	<= 14	15 - 17	>= 18		R S
Ceftiofur	CFU	<= 17	18 - 20	>= 21		R S
Ceftriaxone	CRO	<= 13	14 - 20	>= 21		R S
Cefuroxime = Zinacef®	CXM	<= 14	15 - 17	>= 18		R S
Cephalexin	CFX	<= 14	15 - 17	>= 18		R S
Cephalothin	CF	<= 14	15 - 17	>= 18		R S
Chloramphenicol	C	<= 17	18 - 20	>= 21		R S
Cilastatin&Imipenem = Primaxin®						R S
Ciprofloxacin	CIP	<= 14	16 - 20	>= 21		R S
Clindamycin	CC	<= 13	14 - 17	>= 19		R S
Doxycycline	D	<= 12	13 - 15	>= 16		R S
Enrofloxacin	ENO	<= 17	18 - 20	>= 21		R S
Erythromycin	E	<= 13	14 - 17	>= 18		R S
Gentamicin	GM	<= 6	7 - 9.	>= 10		R S
Imipenem	IPM	<= 13	14 - 15	>= 16		R S
Lincomycin	L	<= 14	15 - 17	>= 18		R S
Marbofloxacin / Zeniquin	MAR	<= 14	15 - 19	>= 20		R S
Meropenem	MEM	<= 13	14 - 15	>= 16		R S
Metronidazole						R S
Neomycin	NE	<= 12	13 - 16	>= 17		R S
Nitrofurantoin	F/M	<= 14	15 - 16	>= 17		R S
Novobiocin	NB	<= 17	18 - 20	>= 21		R S
Ofloxacin	OFX	<= 12	13 - 15	>= 16		R S
Oxacillin	OX	<= 10	11 - 12.	>= 13		R S
Oxytetracycline	T	<= 15	16 - 25	>= 25		R S
Penicillin	P	<= 28	n/a	>= 29		R S
Piperacillin	PIP	<= 17	18 - 20	>= 21		R S
Piperacillin & Tazobactam	TZP	<= 17	18 - 20	>= 21		R S
Pirlimycin HCl = Pirsue®						R S
Polymixin B	PB	<= 8	9 - 11.	>= 12		R S
Rifampin	RA	<= 16	17 - 19	>= 20		R S
SMZ/SXT Tribriessen/Bactrim/Septa	SXT	<= 10	11 - 15.	>= 16		R S
Telavancin	TLV	<= 14	n/a	>= 15		R S
Tetracycline	TE	<= 14	15 - 18	>= 19		R S
Ticarcillin	TC	<= 14	15 - 19	>= 20		R S
Timentin=Ticarcillin&Clavulanic Acid	TIM	<= 15	16 - 20	>= 21		R S
Tobramycin	TN	<= 12	13 - 14	>= 15		R S
Vancomycin	VA	<= 14	15 - 16	>= 17		R S

Report Notes:

To request additional antibiotic, please contact our office. We will work to source the antibiotics that best fit your practice formulary.